

Human IgE Mouse mAb

Cat. No:	V5406
Name:	Human IgE Mouse mAb
Storage:	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	6C3
Host/Source:	Mouse
Subtype:	IgG1
Reactivity:	(predicted: Human)
Purity:	>90% as determined by SDS-PAGE
Format:	Liquid
Buffer:	0.01M PBS (pH7.4).
Applications:	ELISA, ICA
Antigen:	Purified native protein (bs-2226P)
Pairing Information:	

ELISA		Detection		
		V5402	V5404	V5406
Capture	V5402	-	-	-
	V5404	+++	-	-
	V5406	+++	-	-
ICA		Detection		
		V5402	V5404	V5406
Capture	V5402	-	-	-
	V5404	+++	-	-
	V5406	+++	-	-