
Recombinant toxoplasma GRA6 protein

Catalog Number: bs-41556P

Concentration: 1.14mg/ml

Species: Toxoplasma gondii

Predicted MW: 24

Tags: C-His

Purity: >95% as determined by SDS-PAGE

Form: Lyophilized or Liquid

Storage: PBS with 25mM arginine.

Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.

Background: Recombinant toxoplasma GRA6 is a protein with 180 amino acids, used to test the specific IgG and IgM antibody for the diagnosis of Toxoplasma gondii infection, it was expressed from E. coli, containing a 6 x his tag at its C-terminus. This protein usually forms dimer used for ELISA and lateral flow test.