
Recombinant toxoplasma GRA7 protein

Catalog Number: bs-41557P

Concentration: 1.4mg/ml

Species: Toxoplasma gondii

Predicted MW: 26

Tags: GST

Purity: >90% as determined by SDS-PAGE

Form: Lyophilized or Liquid

Storage: 50mM Tris-HCl, pH8.0.

Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.

Background: Recombinant GRA7 is an antigen used to test the specific Toxoplasma gondii antibody, it contains 146 amino acid from GRA7 gene and fused with GST tag.