
Recombinant toxoplasma GRA8 protein

Catalog Number: bs-41558P

Concentration: 0.4mg/ml

Species: *Toxoplasma gondii*

Predicted MW: 20

Tags: C-His

Purity: >95% as determined by SDS-PAGE

Form: Lyophilized or Liquid

Storage: Phosphate buffered saline, 0.05mM EDTA, pH 7.4.

Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.

Background: Recombinant GRA8 is an antigen used to test the specific *Toxoplasma gondii* antibody for the diagnosis of *Toxoplasma gondii* infection. It was expressed from *E. coli*, containing 217 amino acids with 6 x his tag at its C-terminus