
Recombinant *Toxoplasma gondii* SAG2 protein

Catalog Number: bs-41560P

Concentration: 0.56mg/ml

Species: *Toxoplasma gondii*

Predicted MW: 19

Tags: C-His

Purity: >95% as determined by SDS-PAGE

Form: Lyophilized or Liquid

Storage: Phosphate buffered saline.

Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.

Background: *Toxoplasma gondii* is an obligate intracellular protozoan parasite, it has been known associating with severe congenital defects and life threatening infection in immunocompromised hosts. In human, *T. gondii* infection in pregnant women can result in congenital infection, such as late-onset eye disease, it is also a frequent cause of encephalitis in AIDS patients: in animal, the infection can cause a variety diseases and symptoms.