

---

## Recombinant human Endomucin protein, C-His (HEK293)

Catalog Number: bs-43605P

Concentration: >1mg/ml

Species: Human

AA Seq: 19-190/261

Predicted MW: 18.9

Tags: C-His

Endotoxin: Not analyzed

Purity: >90% as determined by SDS-PAGE


Purification: AC

Storage: Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.

Background: Endothelial sialomucin, also called endomucin or mucin-like sialoglycoprotein, which interferes with the assembly of focal adhesion complexes and inhibits interaction between cells and the extracellular matrix.

### VALIDATION IMAGES

---



The purity of the protein is greater than 90% as determined by reducing SDS-PAGE.