

Human LH Mouse mAb

Cat. No:	V6405
Name:	Human LH Mouse mAb
Storage:	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	6H1
Host/Source:	Mouse
Subtype:	IgG
Reactivity:	(predicted: Human)
Format:	Lyophilized or Liquid
Buffer:	0.01M PBS (pH7.4).
Applications:	ELISA
Antigen:	Recombinant Human LH Protein (bs-43083P)
Pairing Information:	

ELISA		Detection					
		V6401	V6402	V6403	V6404	V6405	V6406
Capture	V6401	-	-	-	-	-	-
	V6402	-	-	-	-	-	-
	V6403	-	-	-	+++	-	-
	V6404	-	-	-	-	-	-
	V6405	-	-	-	-	-	+++
	V6406	-	-	-	-	-	-
ICA		Detection					
		V6401	V6402	V6403	V6404	V6405	V6406
Capture	V6401	-	-	-	+++	-	-
	V6402	-	-	-	+	-	-
	V6403	-	-	-	+	-	-
	V6404	-	-	-	-	-	-
	V6405	-	-	-	-	-	-
	V6406	-	-	-	-	-	-