

Human LH Mouse mAb

Cat. No: V6406
Name: Human LH Mouse mAb
Storage: Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type: Monoclonal
Clone No: 18H11
Host/Source: Mouse
Subtype: IgG1
Reactivity: (predicted: Human)
Format: Lyophilized or Liquid
Buffer: 0.01M PBS (pH7.4).
Applications: ELISA
Antigen: Recombinant Human LH Protein ([bs-43083P](#))
Pairing Information:

ELISA		Detection					
		V6401	V6402	V6403	V6404	V6405	V6406
Capture	V6401	-	-	-	-	-	-
	V6402	-	-	-	-	-	-
	V6403	-	-	-	+++	-	-
	V6404	-	-	-	-	-	-
	V6405	-	-	-	-	-	+++
	V6406	-	-	-	-	-	-
ICA		Detection					
		V6401	V6402	V6403	V6404	V6405	V6406
Capture	V6401	-	-	-	+++	-	-
	V6402	-	-	-	+	-	-
	V6403	-	-	-	+	-	-
	V6404	-	-	-	-	-	-
	V6405	-	-	-	-	-	-
	V6406	-	-	-	-	-	-