

TSH Mouse mAb

Cat. No: V6701
 Name: TSH Mouse mAb
 Storage: Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
 Clone type: Monoclonal
 Clone No: 7H1
 Host/Source: Mouse
 Subtype: IgG
 Reactivity: (predicted: Human)
 Format: Lyophilized or Liquid
 Buffer: 0.01M PBS (pH7.4).
 Applications: ELISA
 Antigen: Recombinant Human TSH Protein ([bs-43082P](#))
 Pairing Information:

ELISA		Detection	
		V6701	V6702
Capture	V6701	-	+++
	V6702	-	-