

## TSH Mouse mAb

Cat. No: V6702  
Name: TSH Mouse mAb  
Storage: Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.  
Clone type: Monoclonal  
Clone No: 38B1  
Host/Source: Mouse  
Subtype: IgG1  
Reactivity: (predicted:&nbsp;Human)  
Format: Lyophilized or Liquid  
Buffer: 0.01M PBS (pH7.4).  
Applications: ELISA  
Antigen: Recombinant Human TSH Protein ([bs-43082P](#))  
Pairing Information:

ELISA		Detection	
		V6701	V6702
Capture	V6701	-	+++
	V6702	-	-