

Human IL-6 Mouse mAb

Cat. No:	V0405
Name:	Human IL-6 Mouse mAb
Storage:	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	1E3
Host/Source:	Mouse
Subtype:	IgG2b
Reactivity:	Human
Format:	Lyophilized or Liquid
Buffer:	0.01M PBS (pH7.4).
Applications:	ELISA, ICA
Antigen:	Recombinant Human IL-6 Protein (bs-43079P)
Pairing Information:	

ELISA		Detection				
		V0405	V0406	V0407	V0408	V0409
Capture	V0405	-	+++	-	-	-
	V0406	-	-	-	-	-
	V0407	-	+	-	+++	-
	V0408	-	-	-	-	-
	V0409	-	-	-	-	-
ICA		Detection				
		V0405	V0406	V0407	V0408	V0409
Capture	V0405	-	-	-	-	-
	V0406	-	-	-	-	-
	V0407	+++	-	-	-	+++
	V0408	-	-	-	-	-
	V0409	-	-	-	-	-