

Human IL-6 Mouse mAb

Cat. No: V0407
Name: Human IL-6 Mouse mAb
Storage: Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type: Monoclonal
Clone No: 4C6
Host/Source: Mouse
Subtype: IgG1
Reactivity: Human
Format: Lyophilized or Liquid
Buffer: 0.01M PBS (pH7.4).
Applications: ELISA, ICA
Antigen: Recombinant Human IL-6 Protein ([bs-43079P](#))
Pairing Information:

ELISA		Detection				
		V0405	V0406	V0407	V0408	V0409
Capture	V0405	-	+++	-	-	-
	V0406	-	-	-	-	-
	V0407	-	+	-	+++	-
	V0408	-	-	-	-	-
	V0409	-	-	-	-	-
ICA		Detection				
		V0405	V0406	V0407	V0408	V0409
Capture	V0405	-	-	-	-	-
	V0406	-	-	-	-	-
	V0407	+++	-	-	-	+++
	V0408	-	-	-	-	-
	V0409	-	-	-	-	-