

**bs-0896R**

**[ Primary Antibody ]**

## HHV-5 UL23 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

### — DATASHEET —

**Host:** Rabbit

**Isotype:** IgG

**Clonality:** Polyclonal

**Target:** HHV-5 UL23

**Immunogen:** KLH conjugated synthetic peptide derived from HCMV UL23: 201-284/284.

**Purification:** affinity purified by Protein A

**Concentration:** 1mg/ml

**Storage:** 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.  
Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

**Applications:** **IHC-P** (1:100-500)

**IHC-F** (1:100-500)

**IF** (1:100-500)

**ELISA** (1:5000-10000)

**Reactivity:** (predicted: HHV5)

**Predicted MW.:** 33 kDa