

## Recombinant Mycobacterium bovis MPB70 protein, N-His

Catalog Number: bs-42190P

Concentration: >0.5mg/ml

Species: Mycobacterium bovis

AA Seq: 31-193/193

Predicted MW: 18.6 kDa

Tags: N-His

Endotoxin: Not analyzed

Purity: >90% as determined by SDS-PAGE


Purification: affinity purified by Protein A

Form: Lyophilized or Liquid

Storage: PBS (pH7.4).

Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.

### VALIDATION IMAGES



The purity of the protein is greater than 90% as determined by reducing SDS-PAGE.