
Recombinant human FOXP3 protein, N-His

Catalog Number: bs-42173P

Concentration: >0.5mg/ml

Species: Human

AA Seq: 1-193/431

Predicted MW: 22.2

Detected MW: 25

Tags: N-His

Endotoxin: Not analyzed

Purity: >90% as determined by SDS-PAGE

Purification: affinity purified by Protein A

Form: Lyophilized or Liquid

Storage: 20mM Tris-HCl (pH8.0) with 8M Urea.

Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.

Background: The protein encoded by this gene is a member of the forkhead/winged-helix family of transcriptional regulators. Defects in this gene are the cause of immunodeficiency polyendocrinopathy, enteropathy, X-linked syndrome (IPEX), also known as X-linked autoimmunity-immunodeficiency syndrome. Alternatively spliced transcript variants encoding different isoforms have been identified. [provided by RefSeq, Jul 2008].

VALIDATION IMAGES



The purity of the protein is greater than 80% as determined by reducing SDS-PAGE.