
Human IgE Mouse mAb, Biotin conjugated

Cat. No:	V5405-Bio
Name:	Human IgE Mouse mAb, Biotin conjugated
Storage:	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	1B9
Host/Source:	Mouse
Subtype:	IgG
Reactivity:	(predicted: Human)
Buffer:	PBS (pH=7.4).
Applications:	ELISA
Antigen:	Human IgE