

---

## Human IgE Mouse mAb, Biotin conjugated

Cat. No:	V5406-Bio
Name:	Human IgE Mouse mAb, Biotin conjugated
Storage:	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	6C3-1B7
Host/Source:	Mouse
Subtype:	IgG
Reactivity:	(predicted:&nbsp;Human)
Buffer:	PBS.
Applications:	ELISA, ICA
Antigen:	Human IgE