

www.bioss.com.cn sales@bioss.com.cn techsupport@bioss.com.cn 400-901-9800

zearalenone Mouse mAb

Cat. No: V10801

Name: zearalenone Mouse mAb

Storage: Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Clone type: Monoclonal

Clone No: 4A1

Host/Source: Mouse Subtype: IgG

Reactivity: (predicted: Zearalenone)

Concentration: 1mg/ml

Buffer: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.

Applications: ELISA, ICA

Antigen: Zearalenone conjugated to Bovine Serum Albumin