
vomitoxin (DON) Mouse mAb

Cat. No:	V12201
Name:	vomitoxin (DON) Mouse mAb
Storage:	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	8E5
Host/Source:	Mouse
Subtype:	IgG
Reactivity:	(predicted: vomitoxin)
Concentration:	1mg/ml
Buffer:	0.01M PBS (pH7.4)
Applications:	ELISA, ICA
Antigen:	vomitoxin conjugated to Bovine Serum Albumin