

PGII Mouse mAb

Cat. No: V4201
Name: PGII Mouse mAb
Storage: Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type: Monoclonal
Clone No: 3D5
Host/Source: Mouse
Subtype: IgG2b
Reactivity: Human
Format: Lyophilized or Liquid
Buffer: PBS (pH7.4).
Applications: ELISA
Antigen: Recombinant Human PGII Protein ([bs-43143P](#))
Pairing Information:

ELISA		Detection				
		V4201	V4202	V4203	V4204	V4205
Capture	V4201	-	+++	-	-	-
	V4202	-	-	-	-	-
	V4203	-	+++	-	-	-
	V4204	-	-	-	-	+++
	V4205	-	-	-	-	-