

## PGII Mouse mAb

Cat. No: V4203  
Name: PGII Mouse mAb  
Storage: Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.  
Clone type: Monoclonal  
Clone No: 1B8  
Host/Source: Mouse  
Subtype: IgG2b  
Reactivity: Human  
Format: Lyophilized or Liquid  
Buffer: PBS (pH7.4).  
Applications: ELISA  
Antigen: Recombinant Human PGII Protein ([bs-43143P](#))  
Pairing Information:

| ELISA   |       | Detection |       |       |       |       |
|---------|-------|-----------|-------|-------|-------|-------|
|         |       | V4201     | V4202 | V4203 | V4204 | V4205 |
| Capture | V4201 | -         | +++   | -     | -     | -     |
|         | V4202 | -         | -     | -     | -     | -     |
|         | V4203 | -         | +++   | -     | -     | -     |
|         | V4204 | -         | -     | -     | -     | +++   |
|         | V4205 | -         | -     | -     | -     | -     |