

CA125 Mouse mAb

Cat. No:	V7401
Name:	CA125 Mouse mAb
Storage:	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	10A3
Host/Source:	Mouse
Subtype:	IgG1
Reactivity:	Human
Format:	Lyophilized or Liquid
Buffer:	PBS (pH7.4).
Applications:	ELISA, ICA
Antigen:	Recombinant human MUC16 Protein (bs-43121P)
Pairing Information:	

ELISA		Detection		
		V7401	V7402	V7403
Capture	V7401	-	+++	+++
	V7402	-	-	-
	V7403	-	-	-
ICA		Detection		
		V7401	V7402	V7403
Capture	V7401	-	+++	+++
	V7402	-	-	-
	V7403	-	-	-

PRODUCT SPECIFIC PUBLICATIONS

[IF=34.9] Yi Yakun. et al. Nanopore-based enzyme-linked immunosorbent assay for cancer biomarker detection. NAT NANOTECHNOL. 2025 May;;1-8 ELISA ; Human . 40369343