

CA125 Mouse mAb

Cat. No: V7402
Name: CA125 Mouse mAb
Storage: Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type: Monoclonal
Clone No: 5B12
Host/Source: Mouse
Subtype: IgG2a
Format: Lyophilized or Liquid
Buffer: PBS (pH7.4).
Applications: ELISA, ICA
Antigen: Recombinant human MUC16 Protein ([bs-43121P](#))
Pairing Information:

ELISA		Detection		
		V7401	V7402	V7403
Capture	V7401	-	+++	+++
	V7402	-	-	-
	V7403	-	-	-
ICA		Detection		
		V7401	V7402	V7403
Capture	V7401	-	+++	+++
	V7402	-	-	-
	V7403	-	-	-

PRODUCT SPECIFIC PUBLICATIONS

[IF=34.9] Yi Yakun. et al. Nanopore-based enzyme-linked immunosorbent assay for cancer biomarker detection. NAT NANOTECHNOL. 2025 May;;1-8 ELISA ; Human . 40369343