

## CA125 Mouse mAb

Cat. No: V7403  
Name: CA125 Mouse mAb  
Storage: Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.  
Clone type: Monoclonal  
Clone No: 7A12  
Host/Source: Mouse  
Subtype: IgG2a  
Reactivity: Human  
Format: Lyophilized or Liquid  
Buffer: PBS (pH7.4).  
Applications: ELISA, ICA  
Antigen: Recombinant human MUC16 Protein ([bs-43121P](#))  
Pairing Information:

ELISA		Detection		
		V7401	V7402	V7403
Capture	V7401	-	+++	+++
	V7402	-	-	-
	V7403	-	-	-
ICA		Detection		
		V7401	V7402	V7403
Capture	V7401	-	+++	+++
	V7402	-	-	-
	V7403	-	-	-