
Human IL-6 Mouse mAb, Biotin conjugated

| | |
|---------------|--|
| Cat. No: | V0407-Bio |
| Name: | Human IL-6 Mouse mAb, Biotin conjugated |
| Storage: | Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. |
| Clone type: | Monoclonal |
| Clone No: | 4C6 |
| Host/Source: | Mouse |
| Subtype: | IgG1 |
| Reactivity: | Human |
| Buffer: | 0.01M PBS(pH7.4). |
| Applications: | ELISA, ICA |
| Antigen: | Recombinant Human IL-6 Protein |