

---

## Human IL-6 Mouse mAb, BF488 conjugated

Cat. No:	V0408-BF488
Name:	Human IL-6 Mouse mAb, BF488 conjugated
Storage:	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	4D1
Host/Source:	Mouse
Subtype:	IgG1
Reactivity:	Human
Buffer:	0.01M PBS(pH7.4).
Applications:	ELISA
	Excitation spectrum: 488nm
	Emission spectrum: 519nm
Antigen:	Recombinant Human IL-6 Protein