
beta-Amyloid(1-42) Mouse mAb, PE conjugated

Cat. No:	V5304-PE
Name:	beta-Amyloid(1-42) Mouse mAb, PE conjugated
Storage:	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	3D12
Host/Source:	Mouse
Subtype:	IgG2a
Reactivity:	(predicted: Human)
Buffer:	0.01M PBS(pH7.4).
Applications:	ELISA, ICA Excitation spectrum: 496/564nm Emission spectrum: 578nm
Antigen:	KLH conjugated synthetic peptide derived from human beta-Amyloid(1-42)