

---

## phospho-Tau (Thr217) Mouse mAb, PE conjugated

Cat. No:	V5208-PE
Name:	phospho-Tau (Thr217) Mouse mAb, PE conjugated
Storage:	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	3C4
Host/Source:	Mouse
Subtype:	IgG2b
Reactivity:	(predicted:&nbsp;Human)
Buffer:	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Applications:	ELISA Excitation spectrum: 496nm,564nm Emission spectrum: 578nm Not yet tested in other applications. Optimal working dilutions must be determined by the end user.
Antigen:	KLH conjugated Synthesised phosphopeptide derived from human Tau around the phosphorylation site of Thr217