
phospho-Tau (Thr181) Mouse mAb, PE conjugated

Cat. No:	V5211-PE
Name:	phospho-Tau (Thr181) Mouse mAb, PE conjugated
Storage:	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	12B5
Host/Source:	Mouse
Subtype:	IgG1
Reactivity:	(predicted: Human)
Buffer:	0.01M PBS(pH7.4).
Applications:	ELISA Excitation spectrum: 496/564nm Emission spectrum: 578nm
Antigen:	KLH conjugated Synthesised phosphopeptide derived from human Tau around the phosphorylation site of Thr181