
phospho-Tau (Thr231) Mouse mAb, PE conjugated

Cat. No:	V5223-PE
Name:	phospho-Tau (Thr231) Mouse mAb, PE conjugated
Clone type:	Monoclonal
Clone No:	2D3
Host/Source:	Mouse
Subtype:	IgG1
Reactivity:	(predicted: Human)
Concentration:	1mg/ml
Buffer:	0.01M PBS(pH7.4).
Antigen:	KLH conjugated Synthesised phosphopeptide derived from human Tau around the phosphorylation site of Thr231