

---

## Human Tau Mouse mAb, Biotin conjugated

|               |  |
|---------------|--|
| Cat. No:      | V5213-Bio  |
| Name:         | Human Tau Mouse mAb, Biotin conjugated   |
| Storage:      | Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. |
| Clone type:   | Monoclonal   |
| Clone No:     | 7A2  |
| Host/Source:  | Mouse  |
| Subtype:      | IgG2a  |
| Reactivity:   | (predicted:&nbsp;Human)  |
| Buffer:       | 0.01M PBS(pH7.4).  |
| Applications: | ELISA  |
| Antigen:      | Recombinant human Tau-4 protein  |