

---

## **beta-Amyloid (1-42) Mouse mAb, Biotin conjugated**

Cat. No:	V5306-Bio
Name:	beta-Amyloid (1-42) Mouse mAb, Biotin conjugated
Storage:	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	6H6
Host/Source:	Mouse
Subtype:	IgG2a
Reactivity:	(predicted:&nbsp;Human)
Buffer:	0.01M PBS(pH7.4).
Applications:	ELISA, ICA
Antigen:	KLH conjugated synthetic peptide derived from human beta-Amyloid (1-42)