

---

## Human GFAP Mouse mAb

Cat. No:	V5601
Name:	Human GFAP Mouse mAb
Storage:	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	5D12
Host/Source:	Mouse
Reactivity:	Human
Format:	Liquid
Buffer:	PBS (pH7.4).
Applications:	ELISA
Antigen:	Recombinant Human GFAP Protein ( <a href="#">bs-42001P</a> )
Pairing Information:	
V5601-V5602	