
Recombinant *Proteus hauseri* Endonuclease protein, N-His

Catalog Number: bs-42265P

Concentration: >1mg/ml

Species: *Proteus hauseri*

AA Seq: 1-260/260

Predicted MW: 30.8

Tags: N-His

Endotoxin: Not analyzed

Purity: >90% as determined by SDS-PAGE

Purification: AC

Form: Lyophilized or Liquid

Storage: 20mM Tris-HCl (pH8.0) with 8M Urea.

Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.

VALIDATION IMAGES



The purity of the protein is greater than 90% as determined by reducing SDS-PAGE.