
Prostaglandin F2a-BSA

Catalog Number: bs-2341PB

Form: Lyophilized or Liquid

Storage: 0.01M PBS(pH7.4).

Stored at -20°C. Avoid repeated freeze/thaw cycles.

Background: 生物活性：Dinoprost(Prostaglandin F2_alpha_)是一前列腺素。prostaglandinF2alpha可通过直接作用于黄体素，刺激子宫肌层活性、使宫颈舒张、抑制黄体类固醇生成、诱导黄体溶解。

作用靶点：Prostaglandin Receptor; Endogenous Metabolite; Autophagy; Apoptosis;

作用通路：GPCR/G Protein; Metabolic Enzyme/Protease; Autophagy; Apoptosis;