
ODF3L1 Rabbit pAb

Catalog Number: bs-21133R

Target Protein: ODF3L1

Concentration: 1mg/ml

Form: Liquid

Host: Rabbit

Clonality: Polyclonal

Isotype: IgG

Applications: IHC-P (1:100-500), IHC-F (1:100-500), IF (1:100-500), ICC/IF (1:100-500), ELISA (1:5000-10000)

Reactivity: (predicted:Human, Mouse, Horse)

Predicted MW: 31 kDa

Entrez Gene: 161753

Swiss Prot: Q8IXM7

Source: KLH conjugated synthetic peptide derived from human ODF3L1: 51-150/274.

Purification: affinity purified by Protein A

Storage: Preservative: 0.02% Proclin300, Constituents: 1% BSA, 0.01M PBS, pH7.4.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: The exact function of ODF3L1 remains unknown.