

---

## Human IgA Mouse mAb

Cat. No:	V20103
Name:	Human IgA Mouse mAb
Storage:	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	7H8
Host/Source:	Mouse
Subtype:	IgG1
Reactivity:	Human
Purity:	V20104&V20105-V20103
Format:	Liquid
Buffer:	0.01M PBS(pH7.4)
Applications:	ELISA
Antigen:	Purified native protein ( <a href="#">bs-0474P</a> )