

---

## Human IgM Mouse mAb

Cat. No:	V20302
Name:	Human IgM Mouse mAb
Storage:	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	2B8
Host/Source:	Mouse
Subtype:	IgG2b
Reactivity:	Human
Format:	Liquid
Buffer:	0.01M PBS (pH7.4).
Applications:	ELISA
Antigen:	Purified native protein ( <a href="#">bs-0345P</a> )