
Human IgM Mouse mAb

Cat. No:	V20305
Name:	Human IgM Mouse mAb
Storage:	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	4C4
Host/Source:	Mouse
Subtype:	IgG2a
Reactivity:	Human
Format:	Liquid
Buffer:	0.01M PBS (pH7.4).
Applications:	ELISA
Antigen:	Purified native protein (bs-0345P)