
Human Lambda chain Mouse mAb

| | |
|---------------|---|
| Cat. No: | V20502 |
| Name: | Human Lambda chain Mouse mAb |
| Storage: | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |
| Clone type: | Monoclonal |
| Clone No: | 6F10 |
| Host/Source: | Mouse |
| Subtype: | IgG1 |
| Reactivity: | Human |
| Format: | Liquid |
| Buffer: | 0.01M PBS (pH7.4). |
| Applications: | ELISA |
| Antigen: | Purified native protein (bs-0474P) |