
Human Lambda chain Mouse mAb

Cat. No:	V20504
Name:	Human Lambda chain Mouse mAb
Storage:	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	7D1
Host/Source:	Mouse
Subtype:	IgG1
Reactivity:	Human
Format:	Liquid
Buffer:	0.01M PBS (pH7.4).
Applications:	ELISA
Antigen:	Purified native protein (bs-0345P)