

---

## Human Kappa chain Mouse mAb

Cat. No:	V20401
Name:	Human Kappa chain Mouse mAb
Storage:	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	4C2
Host/Source:	Mouse
Subtype:	IgG2a
Reactivity:	Human
Format:	Liquid
Buffer:	0.01M PBS (pH7.4).
Applications:	ELISA
Antigen:	Purified native protein ( <a href="#">bs-0474P</a> )