

---

## Human Lambda chain Mouse mAb

Cat. No:	V20503
Name:	Human Lambda chain Mouse mAb
Storage:	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	4C6
Host/Source:	Mouse
Subtype:	IgG2b
Reactivity:	Human
Format:	Liquid
Buffer:	0.01M PBS (pH7.4).
Applications:	ELISA
Antigen:	Purified native protein ( <a href="#">bs-0474P</a> )