bs-15438R

[Primary Antibody]

www.bioss.com.cn sales@bioss.com.cn techsupport@bioss.com.cn 400-901-9800

HEATR1 Rabbit pAb

DATASHEET —

Host: Rabbit Isotype: IgG

Clonality: Polyclonal

GenelD: 55127 **SWISS:** Q9H583

Target: HEATR1

Immunogen: KLH conjugated synthetic peptide derived from human HEATR1:

51-150/2144.

Purification: affinity purified by Protein A

Concentration: 1mg/ml

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50%

Glycerol.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated

freeze/thaw cycles.

Applications: IHC-P (1:100-500)

IHC-F (1:100-500) **IF** (1:100-500) ICC/IF (1:100-500) **ELISA** (1:5000-10000)

Reactivity: Mouse (predicted: Human,

Rat, Rabbit, Pig, Sheep, Cow, Dog, Horse)

Predicted MW.: 242 kDa

Subcellular Nucleus

- SELECTED CITATIONS -

• [IF=4.4] Michał Kasperczak, et al. Immunohistological Examination of HEATR1 and SLC27A2 Expression in ccRCC Samples to Evaluate Their Potential as Prognostic Markers—A Preliminary Study. CANCERS. 2025 Jan;17(13):2234 IHC ;Human. 40647532