
RBP4 Rabbit pAb

Cat. No:	W6201
Name:	RBP4 Rabbit pAb
Storage:	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Polyclonal
Clone No:	
Host/Source:	Rabbit
Subtype:	IgG
Reactivity:	Human
Format:	Liquid
Buffer:	0.01M PBS (pH7.4) with 0.02% Proclin300.
Applications:	WB, ELISA, TIA
Antigen:	Recombinant human RBP4 protein (bs-43011P)