
Alpha 1 microglobulin Rabbit pAb

Cat. No:	W6101
Name:	Alpha 1 microglobulin Rabbit pAb
Storage:	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Polyclonal
Clone No:	
Host/Source:	Rabbit
Subtype:	IgG
Format:	Liquid
Buffer:	0.01M PBS (pH7.4) with 0.02% Proclin300.
Applications:	ELISA, TIA
Antigen:	Recombinant human Alpha 1 microglobulin protein (bs-43007P)