

---

## AD7c-NTP Mouse mAb

Cat. No:	V3009
Name:	AD7c-NTP Mouse mAb
Storage:	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	13F6
Host/Source:	Mouse
Reactivity:	Human
Purity:	>90% as determined by SDS-PAGE
Format:	Liquid
Buffer:	0.01M PBS (pH7.4).
Applications:	ELISA
Antigen:	Recombinant Human AD7c-NTP Protein ( <a href="#">bs-41194P</a> )
Pairing Information:	
V3008-V3009	