

**bs-2219R**

**[ Primary Antibody ]**

## **Pectinesterase inhibitor 18 Rabbit pAb**



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

### **— DATASHEET —**

**Host:** Rabbit

**Isotype:** IgG

**Clonality:** Polyclonal

**Target:** Pectinesterase inhibitor 18

**Immunogen:** KLH conjugated synthetic peptide derived from Arabidopsis thaliana (thale cress) pectinesterase inhibitor 18 Pectinesterase: 451-557/557.

**Purification:** affinity purified by Protein A

**Concentration:** 1mg/ml

**Storage:** 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.  
Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

**Applications:** **WB** (1:500-2000)

**IHC-P** (1:100-500)

**IHC-F** (1:100-500)

**IF** (1:100-500)

**ELISA** (1:5000-10000)

**Reactivity:** (predicted: Arabidopsis  
Thaliana)

**Predicted  
MW.:** 58 kDa