

bs-2219R

[Primary Antibody]

Pectinesterase inhibitor 18 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit

Isotype: IgG

Clonality: Polyclonal

Target: Pectinesterase inhibitor 18

Immunogen: KLH conjugated synthetic peptide derived from Arabidopsis thaliana (thale cress) pectinesterase inhibitor 18 Pectinesterase: 451-557/557.

Purification: affinity purified by Protein A

Concentration: 1mg/ml

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Applications: **WB** (1:500-2000)

IHC-P (1:100-500)

IHC-F (1:100-500)

IF (1:100-500)

ELISA (1:5000-10000)

Reactivity: (predicted: Arabidopsis
Thaliana)

**Predicted
MW.:** 58 kDa