

**bs-2249R**

**[ Primary Antibody ]**

## ACA11 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

### — DATASHEET —

**Host:** Rabbit

**Isotype:** IgG

**Clonality:** Polyclonal

**Target:** ACA11

**Immunogen:** KLH conjugated synthetic peptide derived from Arabidopsis thaliana ACA11: 201-300/1025.

**Purification:** affinity purified by Protein A

**Concentration:** 1mg/ml

**Storage:** 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.  
Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

**Applications:** IHC-P (1:100-500)

IHC-F (1:100-500)

IF (1:100-500)

ELISA (1:5000-10000)

**Reactivity:** (predicted: Arabidopsis  
Thaliana)

**Predicted  
MW.:** 112 kDa